

# CATALOG 2021 2022

## ORTHOPLASTY

|                                       |        |
|---------------------------------------|--------|
| <b>SICASIL</b> Poly-addition Silicone | Page 2 |
| <b>REPAIRING</b> Silicone             | Page 3 |
| Finishing <b>VARNISH</b>              | Page 3 |
| <b>APPLICATION</b> Accessories        | Page 3 |
| Orthoplasty <b>ON FABRIC</b>          | Page 3 |

SICASIL POLY-ADDITION SILICONE

**Sicasil Poly-addition Silicone**

Bi-component silicone comes in two doughs of identical consistency:

- Pink (base)
- White (catalyzer).

Mixed together, they make it possible to obtain an orthoplasty product ready to use.

The lot includes two jars of 250 grams each (A+B).

**Instructions**

Pick out the same quantity of product A and B. Mix well for 1 minute. The product is ready to be used.

Polymerization time is about 3 minutes at ambient air temperature of 20 °C (68 °F).

With Hard or Soft Sicasil, grinding can be done just after polymerization. With Super Soft or Extra Soft Sicasil, you have to wait 24 hours after polymerization to grind.



**HARD SICASIL**

Product No.: 0130

**Application**

For all orthoplasty types, especially for correction purposes



**SOFT SICASIL**

Product No.: 0131

**Application**

For protective and comfort orthoplasty



**SUPER SOFT SICASIL**

Product No.: 0132

**Application**

For interdigital orthoplasty, integumentary problems



**EXTRA SOFT SICASIL**

Product No.: 0142

**Application**

For orthoplasty of hyperprotection : perfect for the diabetic or rheumatic feet

**Recommendations**

- Never use the same spoon or the same fingers for both products to avoid polymerization within the jars.
- Never mix the cover A and B.
- Do not wear latex or nitrile gloves.
- You must wear vinyl gloves, and wash carefully your hands with a mild soap to remove powder residue which could interfere with the polymerization.
- Avoid the hydro alcoholic antiseptic solutions.
- Some cosmetic product, such as hydrating cream can harm the polymerization.
- In order to properly store Sicasil, it is important to store it at a temperature below 30 °C.

REPAIRING SILICONE



**SICAFLEX  
REPAIRING  
SILICONE**

Product No.: 0135S



**BATCH OF 5 MIXING  
TIPS**

Product No.: 0135E

FINISHING VARNISH



**VARNISH**

Product No.: 0111S



**BATCH OF 5 MIXING  
TIPS**


Product No.: 0111E

APPLICATION ACCESSORIES



**INJECTION  
SYRINGE**

Product No.: 0139



**APPLICATION TIPS**

Product No.: 16105E

Batch of 100

ORTHOPLASTY ON FABRIC



**ORTHOPLASTY  
WEFT**

Product No.: 0137

Batch of 3 bands -  
50 x 6 cm